

iZoom Loupes

SEEING IS BELIEVING

Affordable high performance loupes provides high resolution and high contrast images for professional use. Custom lightweight frames and optics design deliver comfort and best working posture. Improve efficiency. Reduce fatigue and chronicle back pains.

Part Numbers	DESCRIPTION
8500050-XXX_XC	Adilte Frame
8500051-XXX_XC	Metal Frame
8500052-YYYY_XC	Variable Rudilite Frame

Including soft case, cleaning cloth, lens cover, headstrap and screw driver.

Magnification XXX: 200 - 2.0X; 250 - 2.5X; 300 - 3.0X; 400 - 4X;
 Magnification YYYY: 2530 - 2.5 to 3.0; 3035 - 3.0 to 3.5 variable magnification
 Working Distance: _X: S (short), R (regular), L(long), XL (extra long)

Working Distance(inch)	S	R	L	XL
	11-14	14-18	17-21	20-24



HD TTL loupes are economic lightweight 2.5x magnification device compatible with Nano light for the most comfortable operation.

Working Distance: 14-20 inches
 Field of view: 4 inches
 Depth of field: 5-6 inches

Part Numbers	DESCRIPTION
8500054-XXS	2.5X Custom TTL (Rudilite frame)

Metal Case Included; Prescription insert optional



DentLight Inc.
 1825 Summit Ave. #210
 Plano, TX 75074, USA
 sales@dentlight.com
 www.dentlight.com