



Natural GP™

The Natural GP™ gutta percha used for the Soft-Core® Endodontic Obturator is processed from natural raw rubber harvested from Indonesian rubber trees.

All manufacturing is performed by CMS Dental A/S in Denmark to exact specifications.

The ingredients of the Natural GP™ are free of any toxic components.

Each step of manufacturing for every batch of Natural GP™ is stringently controlled to ensure optimal product quality.

Soft-Core® Endodontic Obturator

– as good as it gets!



Agriculture rubber plantation early in the morning

Detachable handle

New improved heater

Compatible with all file systems

More than 25 million sold worldwide

Quality since 1996

Soft-Core® Endodontic Obturator

A radiopaque polymere core with a thermoplastic gutta percha

- Flexibility makes even the most challenging root canals possible to fill
- Quick, efficient, reliable and filled in seconds
- Compatible with any rotary file system
- A compatible match with any tapered preparation of the root canal
- Has special features for filling in narrow, curved and long canals
- Easy to use and convenient snap off excess gutta percha
- Clinically proven for more than 18 years
- Safe

Soft-Core®
Endodontic
Obturator #25

Available in sizes:
#20, #25, #30, #35,
#40, #45, #50, #55,
#60 and #70.



Soft-Core® Heater



Removing the foam
before leaving the
mixture to dry

Did you know....

Soft-Core®'s carrier is **thinner** than other obturator systems, allowing greater flexibility to navigate curved canals.

A thinner core also means **more gutta percha** for denser fills every time. And with just **one universal taper size**, you can greatly reduce and simplify your carrier inventory.

The efficient choice for a predictable 3-dimensional filling

Over the years we tried to improve the Soft-Core® Endodontic Obturators but again and again we came to the conclusion - *it's as good as it gets!*

Yet we decided to improve the Soft-Core® Heater and achieved:

- Instant heating
- Heats 4 obturators at the same time
- Keeps obturators warm and ready for use for 15 minutes
- Root filling mission accomplished in less than half the time, usually spent filling a root canal

A 3D filled root canal means that lateral and accessory canals are filled. The unique design of the Soft-Core® Endodontic Obturator makes it the preferred choice for several dentists.

It all begins with the patented design of the carrier with a detachable handle - the handle can be bend to ease treatment of posterior teeth.

The carrier is covered with flowable gutta percha. One single obturator will fill the entire canal and because of the hydraulic force, the warm flowable gutta percha is carried to the apex of the tooth.

The CMS Dental obturator has a slight oval carrier which provides an excellent backflow of the gutta percha coronally.