





endodontics

CHANGZHOU BOMEDENT MEDICAL TECHNOLOGY CO.,LTD
Add:No.9 Changyang Road West Taihu Science &
Technology industrial park Changzhou,Jiangsu,China
Email: wu.feib@bomedent.com
www.bomedent.com
Tel: 0086-0519-89881606

CHANGZHOU BOMEDENT MEDICAL TECHNOLOGY CO.,LTD is a high-tech enterprise who focus on researching, developing, selling and service dental instruments. Our main products are endodontics micromotor, apex locator and related products. we will work closely with purpose to establish a new domestic oral medical environment, research and develop, produce, sales high-end oral medical equipment, we own a R & D team who has more than ten years experience in the development of dental products. Our team will take the technology as the core, depend on the quality as life, keep continuous innovation, dedicate to provide the highest cost-effective products and good after service.



we are located in the beautiful city ChangZhou, Yangtze River delta



<https://stomshop.pro/>

ENDO MOTOR



Micro-motor for root canal preparation with intergrated apex locator


iRoot pro

➤ **APEX LOCATOR** + **ENDO MOTOR**



pink color for option


FOUR WORKING MODES



type	portable
memory	10 memory setting
rotation	clockwise/anti-clockwise/reciprocating
torque control	auto reverse/auto stop
battery	rechargeable/replaceable
charging	charging base/Micro USB
display	OLED/Blue Tooth
available contra angle	16:1
apex locator control motor working mode	auto start/auto stop/auto reverse

The most convenient root canal treatment system , complete root canal length measurement and root canal preparation on the handpiece, you just need to connect a lip hook with the motor to make the canal length measurement and preparation at the same time

Apex locator control motor working mode: Auto start/Auto stop/Auto reverse



We are committed to producing safe and reliable medical equipment to let doctors use more handy, so we recommend doctors to use the apex locator control motor working mode



The most convenient root canal treatment system, complete root canal length measurement and root canal preparation on the handpiece, you just need to connect a lip hook with the motor to make the canal length measurement and preparation at the same time







Early generation motor requires a tip hook, sleeve and cable complicated and stumble



BOMEDENT motor, no accessories required, clean and clear



Portable(rechargeable), smart, cordless




6 angles available both hands are applicable

INTEGRATED BLUETOOTH 



Bluetooth function, displayed data synchronize with our Endo Motor iRoot pro



Compatible with most NI-TI file on market, flexible



Push button, convenient



High temperature autoclavable parts

iRoot apex



SMART

PRECISE



High resolution, wide angle of view LCD, the root canal changes can be clearly observed from different angles
Based on multiple frequency network impedance measurement technology, automatic calibrating ensures the measurements are accurate.



touch panel



Equipped with high precision gyroscope, switch the display mode according to the state of display



Multi angle adjustment




Optional Bluetooth function, displayed data synchronize with our apex locator iRoot apex, perfect visibility



MARC III



Both hands are applicable



Marc III can work together with iRoot Apex by cable



type	portable
memory	10 memory setting
rotation	clockwise/anti-clockwise/reciprocating
torque control	auto reverse/auto stop
battery	rechargeable/replaceable
charging	charging base/Micro USB
display	OLED
available contra angle	16:1



Push button, convenient



Compatible with most NI-TI file on market , flexible