

Best Choice for Best Treatment!

Dr's Light AT is a new conceptual curing light applied many advantages for user's demands.
Press the button and experience new light technology with Dr's Light AT.



- IC-WHT60
- IC-WHW61
- IC-DCAM70
- IC-DCAU71
- IC-WHC50
- IC-EZAU81

- CL-DLT30
- CL-DL10
- CL-DC20
- CL-AT24

- AL-DF10

'Best Choice for Best Treatment'

Dr's Light AT is a new conceptual curing light applied many advantages for user's demands.
Press the button and experience new light technology with Dr's Light AT.

Features

- Focused light
- 4 curing light modes
- Time setting available (Soft Start is fixed for 20 sec.)
- 330° head twist
- Interchangeable battery pack
- Big capacity battery



Without focus lens



With focus lens



Slim head



Disposable cover





Interchangeable battery


Specification

Model	CL-AT24
Dimensions (Included Battery)	Approx. 232 X 26 X 26 mm
Weight	Approx. 130 g
Wave length	440~490 (nm)
Power source	Lithium-ion Battery
Light source	High power LED

Accessory



Comparison with Other products



Various curing modes

Selectable light intensity and time as you want

