

Dr's Cam

Enjoy Wireless & Auto Focus!

Wireless Autofocus Intraoral Camera with Mouse



Wired type

Enjoy Wireless & Auto Focus!

Wireless Autofocus Intraoral Camera with Mouse

1 Features of Dr's cam

Wireless type only

-  **High speed digital wireless**
High speed digital wireless provide more clear live image without interruption.
-  **Wireless receiver**
Compact wireless receiver without any power supply.
-  **Macro to Full arch**
Camera is able to focus from macro mode (close) to full arch.
-  **Auto focus**
Can get clear image with auto focus function in oral cavity.
-  **HD image**
Can get high definition image without distortion.
-  **Shadowless Lighting**
Rounded 12pcs LEDs make shadowless lighting.
-  **Gyro mouse built-in**
Useable as a mouse.
-  **280° twist head**
Easy to get approach anywhere in oral cavity.
-  **OLED screen**
Bright and clear OLED screen indicates various device status.
-  **Adjustable LED brightness**
To avoid image quality loss by light reflection, can adjust LED brightness.
-  **3 Color Impressions**
Can change live image color impression (Original, Warm, Cool)
-  **Software compatibility**
Can be use camera with most of dental imaging software.
* It would need twain driver for some of dental imaging software.




2 Specification

Model name	IC-DCAM70(Wireless model)	
Camera	Operation voltage	3.7V
	Size	Approx. 255 x 30 x 28 mm
	Weight	Approx. 89g
	Resolution	HD(720P)
	Sensor/Lens type	CMOS / Liquid lens
Charger	Operation voltage	5V
	Size	Approx. 110 x 105 x 68 mm
	Weight	Approx. 100g
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

Model name	IC-DCAU71(Wired model)	
Camera	Operation voltage	5V
	Size	Approx. 202 x 30 x 28 mm
	Weight	Approx. 56g
	Resolution	HD(720P)
	Sensor/Lens type	CMOS / Liquid lens
USB cable	Length	Approx. 3M
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

3 Minimum System Requirement

PC system:	2.0 Ghz CPU, 1GB RAM
OS:	Windows XP service pack 3
USB port:	USB 2.0



Upper images are take by real product.

4 DrsViewer2

Image Viewer included.

