

Модельный ряд шумопоглощающих кожухов «ROMAX».

Таблица-помощник в подборе для Вашего компрессора.

модель компрессора		Romax-K0131	Romax-K0132	Romax-K0133	Romax-K0134
Габариты кожуха, мм		600 x 600 x 850	850 x 700 x 950	1000 x 650 x 900	1300 x 750 x 1050
Вентилятор, шт.		1	2	2	3
Розетка, шт.		1	1	нет	нет
Люк		нет	1	1	1
Мощность, Вт		12	25	25	38
Еком (Словакия)	DK50-10 Z	✓			
	DK50-10 Z/M	✓			
	DK50 PLUS		✓		
	DK50 PLUS/M		✓		
	DK50 2V		✓		
	DK50 2V/M		✓		
	DK50 2V/50		✓		
	DK50 2V/50/M		✓		
	DK50 4VR/50		✓		
	DK50 4VR/50/M		✓		
	DK50 2x2V/110				✓
	DK50 2x2V/110/M				✓
	DK50 2x4VR/110				✓
Cattani (Италия)	Cattani 30-67		✓		
	Cattani 30-160		✓		
	Cattani 45-165		✓		
	Cattani 45-238		✓		
	Cattani 100-320				✓
Remeza (Белоруссия)	KM-24.OLD10	✓			
	KM-24.OLD10Д		✓		
	KM-24.OLD15	✓			
	KM-24.OLD15Д	✓			
	KM-24.OLD20	✓			
	KM-24.GMS 150	✓			
	KM-24.GMS 150Д		✓		
	KM-24.VS 204	✓			
	KM-24.VS 204Д		✓		
	KM-50.VS 204		✓		
	KM-50.VS 204Д		✓		
	KM-50.OLD 15		✓		
	KM-50.OLD 15Д		✓		
	KM-50.OLD 20		✓		
	KM-50.OLD 20Д		✓		
	KM-100.VS 204Т				✓
	KM-100.VS 204ТД				✓
	KM-100.OLD 15Т				✓
	KM-100.OLD 15ТД				✓
	KM-100.OLD 20Т				✓
	KM-100.OLD 20ТД				✓
Mercury (Китай)	HK-1EW-30	✓			
	HK-2EW-35		✓		
	HK-3EW-55			✓	
	HK-4EW-65			✓	
	HK-6EW-90				✓
Eluan (Китай)	JYK35	✓			
	JYK50			✓	
	JYK50A			✓	
	JYK65			✓	
	JYK65A			✓	
	JYK135				✓
	JWA50A			✓	
	JWA30			✓	
	EYK135				✓
	EYK50			✓	
	EYK65			✓	
	EYK35			✓	
	JWA50			✓	
Dental League (Китай)	DL300/2EW-35L		✓		