


Ultrasurgery Tip easy to use instructions

PSK-1 (US1)

MODE: BONE


High effectiveness osteotomy of large bone sections during maxillofacial surgery, it can also be used for exodontia exodontics.



PSK-2 (US2)

MODE: BONE


High effectiveness osteotomy technique in maxilla and mandible (ridge expansion, corticotomy technique, bone block grafting).



PSK-4 (US4)

MODE: BONE


Universal osteoplasty: periodontal osteotomy, bone chips harvesting, inflammatory tissue removal (cyst, etc.).



PSK-5 (US5)

MODE: BONE


High efficiency bone osteoplasty: bone remodelling and harvesting of bone chips



PCK-1 (UC1)

MODE: BONE

Used to cut off the ankylosis and root fraction techniques.



PLK-3 (UL3)

MODE: BONE

Diamond-coated(100 μ m) instrument for micrometric osteotomy or osteoplasty: non-traumatic, to finalize the osteotomy or osteoplasty on thin bone and/or near delicate anatomic structures.

