



Shenzhen Yurucheng Dental Materials Co.,ltd

1 Floor,Building 4,No.1 Branch of Zhongcheng Life Science Park,No.14
Zhongxing Road,Kengzi Street,Pingshan District,Shenzhen, China

Web : www.yucera.com

Tel : 400-995-8505

Email : info@yucera.com

Fax : +86(755)-28228526



Zirconia material
processing solution


5 Axis Dental Milling Machine

TECHNOLOGY INNOVATION
PEOPLE ORIENTED
CONCENTRATING ON R&D



PRODUCT DESCRIPTION >>>

Yucera dental milling machine is a dry cutting equipment with five-axis linkage, designed with micron-level high-precision cutting and high-reliability performance to provide you with a more efficient, powerful, accurate and stable processing experience.



TECHNICAL PARAMETER

Equipment size	530*650*750mm	Processing methods	Five-axis linkage, dry milling
Weight	102kg	Burs specifications	Special props for handle diameter 4mm
Rotating axis operating angle	A: 360° B: ±30°	Main axis speed	0-60,000rpm
Total power	1.5kw	Working pressure	4.5-7.5bar (no water, no gasoline)
Cutting precision	0.2mm	Stable voltage	220-230V
Tool magazine capacity	7	Stable air pressure	≥6.0bar

MACHINABLE MATERIAL >>>



Zirconia Block



PMMA




Wax



Composite Material

PRODUCT STRUCTURE >>>



MACHINABLE PRODUCTS >>>



Veneer Full crown Full crown bridge Inner crown Inner crown bridge



Inlay Telescopic crown Orthodontic braces Complete denture Dental mold

* The above are just a few main products, for more information, please contact us.